



MADIHA ORTHO

Solution for humanity



"Your muscles cannot get 'longer' without some rather radical orthopedic surgery."



CATALOGUE

2021

SMALL FRAGMENT
LOCKING SCREWS

CONTENT

S. No.	Products Name	Page No
1	Locking Screw ø 2.4mm - Variable Angle (Star Like)	3
2	Cortex Screw ø 2.4mm (Star Like) - Self Tapping	3
3	Cortex Screw ø 2.7mm (Star Like) - Self Tapping	4
4	Locking Screw ø 3.5mm - Self Tapping	4
5	Locking Screw ø 3.5mm Self Tapping and Self Drilling	6
6	Locking Screw ø 2.7mm - Self Tapping	7
7	Locking Screw ø 2.4mm - Self Tapping	8



Locking Implants - Small Fragment Screws


Screws

Locking Screw Ø 2.4mm - Variable Angle (Star Like)




Codes		
Length in mm	Stainless Steel	Titanium
12	MOPL0319-01	Ti-MOPL0319-01
14	MOPL0319-02	Ti-MOPL0319-02
16	MOPL0319-03	Ti-MOPL0319-03
18	MOPL0319-04	Ti-MOPL0319-04
20	MOPL0319-05	Ti-MOPL0319-05
22	MOPL0319-06	Ti-MOPL0319-06
24	MOPL0319-07	Ti-MOPL0319-07
26	MOPL0319-08	Ti-MOPL0319-08
82	MOPL0319-09	Ti-MOPL0319-09
84	MOPL0319-10	Ti-MOPL0319-10

Cortex Screw Ø 2.4mm (Star Like) - Self Tapping



Codes		
Length in mm	Stainless Steel	Titanium
12	MOPL0320-01	Ti-MOPL0320-01
14	MOPL0320-02	Ti-MOPL0320-02
16	MOPL0320-03	Ti-MOPL0320-03
18	MOPL0320-04	Ti-MOPL0320-04
20	MOPL0320-05	Ti-MOPL0320-05
22	MOPL0320-06	Ti-MOPL0320-06
24	MOPL0320-07	Ti-MOPL0320-07
26	MOPL0320-08	Ti-MOPL0320-08
20	MOPL0320-09	Ti-MOPL0320-09
84	MOPL0320-10	Ti-MOPL0320-10


Cortex Screw Ø 2.7mm (Star Like) - Self Tapping



Codes		
Length in mm	Stainless Steel	Titanium
12	MOPL0321-01	Ti-MOPL0321-01
14	MOPL0321-02	Ti-MOPL0321-02
16	MOPL0321-03	Ti-MOPL0321-03
18	MOPL0321-04	Ti-MOPL0321-04
20	MOPL0321-05	Ti-MOPL0321-05
22	MOPL0321-06	Ti-MOPL0321-06
24	MOPL0321-07	Ti-MOPL0321-07
26	MOPL0321-08	Ti-MOPL0321-08
82	MOPL0321-09	Ti-MOPL0321-09
84	MOPL0321-10	Ti-MOPL0321-10

Locking Screw Ø 3.5mm - Self Tapping

- Locking Screw Ø 3.5mm- Self Tapping is used for the fixation of Small Fragment Locking Plate.
- Larger core dia of the screw provides strength and threaded head for locking in plate.
- Self tapping eliminates the use of bone tap.
- Drill Bit: 2.8mm




Codes		
Length in mm	Stainless Steel	Titanium
10	MOPL0322-01	Ti-MOPL0322-01
12	MOPL0322-02	Ti-MOPL0322-02

14	MOPL0322-03	Ti-MOPL0322-03
16	MOPL0322-04	Ti-MOPL0322-04
18	MOPL0322-05	Ti-MOPL0322-05
20	MOPL0322-06	Ti-MOPL0322-06
22	MOPL0322-07	Ti-MOPL0322-07
24	MOPL0322-08	Ti-MOPL0322-08
26	MOPL0322-09	Ti-MOPL0322-09
82	MOPL0322-10	Ti-MOPL0322-10
22	MOPL0322-11	Ti-MOPL0322-11
32	MOPL0322-12	Ti-MOPL0322-12
34	MOPL0322-13	Ti-MOPL0322-13
36	MOPL0322-14	Ti-MOPL0322-14
38	MOPL0322-15	Ti-MOPL0322-15
40	MOPL0322-16	Ti-MOPL0322-16
42	MOPL0322-17	Ti-MOPL0322-17
44	MOPL0322-18	Ti-MOPL0322-18
46	MOPL0322-19	Ti-MOPL0322-19
48	MOPL0322-20	Ti-MOPL0322-20
50	MOPL0322-21	Ti-MOPL0322-21
52	MOPL0322-22	Ti-MOPL0322-22
54	MOPL0322-23	Ti-MOPL0322-23
56	MOPL0322-24	Ti-MOPL0322-24
58	MOPL0322-25	Ti-MOPL0322-25
60	MOPL0322-26	Ti-MOPL0322-26

Locking Screw Ø 3.5mm Self Tapping and Self Drilling

→ Drill Bit: 2.8mm




Codes		
Length in mm	Stainless Steel	Titanium
10	MOPL0323-01	Ti-MOPL0323-01
12	MOPL0323-02	Ti-MOPL0323-02
14	MOPL0323-03	Ti-MOPL0323-03
16	MOPL0323-04	Ti-MOPL0323-04
18	MOPL0323-05	Ti-MOPL0323-05
20	MOPL0323-06	Ti-MOPL0323-06
22	MOPL0323-07	Ti-MOPL0323-07
24	MOPL0323-08	Ti-MOPL0323-08
26	MOPL0323-09	Ti-MOPL0323-09
28	MOPL0323-10	Ti-MOPL0323-10
30	MOPL0323-11	Ti-MOPL0323-11
32	MOPL0323-12	Ti-MOPL0323-12
34	MOPL0323-13	Ti-MOPL0323-13
36	MOPL0323-14	Ti-MOPL0323-14
38	MOPL0323-15	Ti-MOPL0323-15
40	MOPL0323-16	Ti-MOPL0323-16
42	MOPL0323-17	Ti-MOPL0323-17
44	MOPL0323-18	Ti-MOPL0323-18
46	MOPL0323-19	Ti-MOPL0323-19
48	MOPL0323-20	Ti-MOPL0323-20

50	MOPL0323-21	Ti-MOPL0323-21
52	MOPL0323-22	Ti-MOPL0323-22
54	MOPL0323-23	Ti-MOPL0323-23
56	MOPL0323-24	Ti-MOPL0323-24
58	MOPL0323-25	Ti-MOPL0323-25
60	MOPL0323-26	Ti-MOPL0323-26

Locking Screw Ø 2.7mm - Self Tapping

- Locking Screw Ø 2.7mm - Self Tapping is used for the Fixation of Small Fragment Locking Plate.
- Larger core dia of the screw provides strength and threaded head for locking in plate.
- Self tapping feature eliminates the use of bone tap.
- Drill Bit: 2.0mm




Code		
Length in mm	Stainless Steel	Titanium
20	MOPL0324-01	Ti-MOPL0324-01
22	MOPL0324-02	Ti-MOPL0324-02
24	MOPL0324-03	Ti-MOPL0324-03
26	MOPL0324-04	Ti-MOPL0324-04
82	MOPL0324-05	Ti-MOPL0324-05
84	MOPL0324-06	Ti-MOPL0324-06
32	MOPL0324-07	Ti-MOPL0324-07
34	MOPL0324-08	Ti-MOPL0324-08
36	MOPL0324-09	Ti-MOPL0324-09
38	MOPL0324-10	Ti-MOPL0324-10
40	MOPL0324-11	Ti-MOPL0324-11
42	MOPL0324-12	Ti-MOPL0324-12
44	MOPL0324-13	Ti-MOPL0324-13
46	MOPL0324-14	Ti-MOPL0324-14

48	MOPL0324-15	Ti-MOPL0324-15
50	MOPL0324-16	Ti-MOPL0324-16
52	MOPL0324-17	Ti-MOPL0324-17
54	MOPL0324-18	Ti-MOPL0324-18
56	MOPL0324-19	Ti-MOPL0324-19
58	MOPL0324-20	Ti-MOPL0324-20
60	MOPL0324-21	Ti-MOPL0324-21

Locking Screw Ø 2.4mm - Self Tapping

→ Drill Bit: 1.8mm



Codes		
Length in mm	Stainless Steel	Titanium
6	MOPL0325-01	Ti-MOPL0325-01
8	MOPL0325-02	Ti-MOPL0325-02
10	MOPL0325-03	Ti-MOPL0325-03
12	MOPL0325-04	Ti-MOPL0325-04
14	MOPL0325-05	Ti-MOPL0325-05
16	MOPL0325-06	Ti-MOPL0325-06
18	MOPL0325-07	Ti-MOPL0325-07
20	MOPL0325-08	Ti-MOPL0325-08
22	MOPL0325-09	Ti-MOPL0325-09
24	MOPL0325-10	Ti-MOPL0325-10
26	MOPL0325-11	Ti-MOPL0325-11
82	MOPL0325-12	Ti-MOPL0325-12
84	MOPL0325-13	Ti-MOPL0325-13



MADIHA ORTHO

Solution for humanity



support@madahaortho.com

info@madahaortho.com

sales@madahaortho.com

Contact: +918448115377
+917008406585

WhatsApp: +918697716553
+917602449449

www.madihaortho.com

Office: 542,4th floor, G.T. Road, Howrah-71101
(West Bengal, India)