




MADIHA ORTHO

Solution for humanity



"Your muscles
cannot get
"longer"
without some
rather radical
orthopedic
surgery."



CATALOGUE

2021

CANNULATED SCREWS

CONTENT

S. No.	Products Name	Page No
1	Cannulated Screw Ø 2.7mm Short Thread	3
2	Washer - S.S.	4
3	Cannulated Screw Ø 2.7mm Full Thread	4
4	Washer - Titanium	5
5	Cannulated Cancellous Screw Ø 3.0mm Short Thread	5
6	Cannulated Cancellous Screw Ø 3.0mm Full Thread	7
7	Cannulated Screw Ø 3.0mm - Short Thread	8
8	Cannulated Screw Ø 3.0mm - Full Thread	10
9	Cannulated Screw Ø 3.5mm Short Thread	12
10	Cannulated Screw Ø 3.5mm Full Thread	13
11	Cannulated Cancellous Screw Ø 4.0mm - Short Thread	14
12	Cannulated Cancellous Screw Ø 4.0mm - Full Thread	15
13	Cannulated Screw Ø 4.5mm Short Thread	17
14	Cannulated Screw Ø 4.5mm - Full Thread	18
15	Cancellous Cannulated Screw Ø 6.5mm 16mm Thread	20
16	Cancellous Cannulated Screw Ø 6.5mm 32mm Thread	21
17	Cancellous Cannulated Screw Ø 6.5mm Full Thread	23
18	Cancellous Cannulated Screw Ø 7.0mm, 16mm Thread	24
19	Cancellous Cannulated Screw Ø 7.0mm, 32mm Thread	25
20	Cancellous Cannulated Screw Ø 7.0mm, Full Thread	26
21	Cancellous Cannulated Screw Ø 7.3mm - 16mm Thread	27
22	Cancellous Cannulated Screw Ø 7.3mm - 32mm Thread	28
23	Cancellous Cannulated Screw Ø 7.3mm - Full Thread	29

CANNULATED SCREWS


Small Cannulated Screws (Ø3.0, Ø3.5, Ø4.0mm) are intended for fixation of intraarticular and extra-articular fractures and on-unions of small bones and small bone fragments; arthrodeses of small joints; bunionectomies and osteotomies; scaphoid and other carpal bones, metacarpals, tarsals, metatarsals, patella, ulnar styloid, capitellum, radial head and radial styloid. Large Cannulated Screws (Ø6.5, Ø7.0, Ø7.3mm) are indicated for the fixation of small and long bone fractures in both adult and pediatric patients, including fractures of the femur, tibia, pelvis, acetabulum, the ulna, bones of the mid foot and the calcaneus.

SCREW

2.7mm


Cannulated Screw Ø 2.7mm Short Thread

→ Drill Bit : 2.0mm



Length in mm	Code - S.S	Code - TITANIUM
	Self Drilling	Self Drilling
17	MOPL0377-01	Ti-MOPL0377-01
18	MOPL0377-02	Ti-MOPL0377-02
19	MOPL0377-03	Ti-MOPL0377-03
20	MOPL0377-04	Ti-MOPL0377-04
21	MOPL0377-05	Ti-MOPL0377-05
22	MOPL0377-06	Ti-MOPL0377-06
23	MOPL0377-07	Ti-MOPL0377-07
24	MOPL0377-08	Ti-MOPL0377-08
25	MOPL0377-09	Ti-MOPL0377-09
26	MOPL0377-10	Ti-MOPL0377-10
27	MOPL0377-11	Ti-MOPL0377-11
28	MOPL0377-12	Ti-MOPL0377-12
29	MOPL0377-13	Ti-MOPL0377-13
30	MOPL0377-14	Ti-MOPL0377-14


Washer - S.S.



Code	Description
MOPL0378-01	Washer for 4.5mm Cannulated Screws, SS 316L
MOPL0378-02	Washer for 6.5mm Cannulated Screws, SS 316L
MOPL0378-03	Washer for 7.0mm Cancellous Screws, SS 316L
MOPL0378-04	Washer for 7.0mm Cannulated Screws, SS 316L
MOPL0378-05	Washer for 7.3mm Screws, SS 316L
MOPL0378-06	Screw Washer for 4.5mm Cannulated Screws, Dia :15mm, SS 316L

Cannulated Screw Ø 2.7mm Full Thread

→ Drill Bit : 2.0mm



Code - S.S		Code - TITANIUM
Length in mm	Self Drilling	Self Drilling
17	MOPL0379-01	Ti-MOPL0379-01
18	MOPL0379-02	Ti-MOPL0379-02
19	MOPL0379-03	Ti-MOPL0379-03
20	MOPL0379-04	Ti-MOPL0379-04
21	MOPL0379-05	Ti-MOPL0379-05
22	MOPL0379-06	Ti-MOPL0379-06
23	MOPL0379-07	Ti-MOPL0379-07
24	MOPL0379-08	Ti-MOPL0379-08
25	MOPL0379-09	Ti-MOPL0379-09
26	MOPL0379-10	Ti-MOPL0379-10
27	MOPL0379-11	Ti-MOPL0379-11
28	MOPL0379-12	Ti-MOPL0379-12
29	MOPL0379-13	Ti-MOPL0379-13
30	MOPL0379-14	Ti-MOPL0379-14

Washer - Titanium




Code - Description

Ti-MOPL0380-01	- Washer for 6.5mm Cannulated Screws, Titanium
Ti-MOPL0380-02	- Washer for 4.5mm Cannulated Screws, Titanium
Ti-MOPL0380-03	- Washer for 7.0mm Cancellous Screws, Titanium
Ti-MOPL0380-04	- Washer for 7.0mm Cannulated Screws, Titanium
Ti-MOPL0380-05	- Washer for 7.3mm Screws, Titanium
Ti-MOPL0380-06	- Screw Washer for 4.5mm Cannulated Screws, Dia :15mm, Titanium


3.0mm

Cannulated Cancellous Screw Ø 3.0mm Short Thread


→ Drill Bit : 2.0mm



Self Drilling
Tip - Titanium



Self Drilling Tip - s.s.




Code - S.S		Code - TITANIUM
Length in mm	Self Drilling	Self Drilling
8	MOPL0381-01	Ti-MOPL0381-01
9	MOPL0381-02	Ti-MOPL0381-02
10	MOPL0381-03	Ti-MOPL0381-03
11	MOPL0381-04	Ti-MOPL0381-04

12	MOPL0381-05	Ti-MOPL0381-05
13	MOPL0381-06	Ti-MOPL0381-06
14	MOPL0381-07	Ti-MOPL0381-07
15	MOPL0381-08	Ti-MOPL0381-08
16	MOPL0381-09	Ti-MOPL0381-09
17	MOPL0381-10	Ti-MOPL0381-10
18	MOPL0381-11	Ti-MOPL0381-11
19	MOPL0381-12	Ti-MOPL0381-12
20	MOPL0381-13	Ti-MOPL0381-13
21	MOPL0381-14	Ti-MOPL0381-14
22	MOPL0381-15	Ti-MOPL0381-15
23	MOPL0381-16	Ti-MOPL0381-16
24	MOPL0381-17	Ti-MOPL0381-17
25	MOPL0381-18	Ti-MOPL0381-18
26	MOPL0381-19	Ti-MOPL0381-19
27	MOPL0381-20	Ti-MOPL0381-20
28	MOPL0381-21	Ti-MOPL0381-21
29	MOPL0381-22	Ti-MOPL0381-22
30	MOPL0381-23	Ti-MOPL0381-23
32	MOPL0381-24	Ti-MOPL0381-24
34	MOPL0381-25	Ti-MOPL0381-25
36	MOPL0381-26	Ti-MOPL0381-26
38	MOPL0381-27	Ti-MOPL0381-27
40	MOPL0381-28	Ti-MOPL0381-28
42	MOPL0381-29	Ti-MOPL0381-29
44	MOPL0381-30	Ti-MOPL0381-30
46	MOPL0381-31	Ti-MOPL0381-31
48	MOPL0381-32	Ti-MOPL0381-32
50	MOPL0381-33	Ti-MOPL0381-33

Cannulated Cancellous Screw Ø 3.0mm Full Thread

→ Drill Bit : 2.0mm




Code - S.S		Code - TITANIUM
Length in mm	Self Drilling	Self Drilling
8	MOPL0382-01	Ti-MOPL0382-01
9	MOPL0382-02	Ti-MOPL0382-02
10	MOPL0382-03	Ti-MOPL0382-03
11	MOPL0382-04	Ti-MOPL0382-04
12	MOPL0382-05	Ti-MOPL0382-05
13	MOPL0382-06	Ti-MOPL0382-06
14	MOPL0382-07	Ti-MOPL0382-07
15	MOPL0382-08	Ti-MOPL0382-08
16	MOPL0382-09	Ti-MOPL0382-09
17	MOPL0382-10	Ti-MOPL0382-10
18	MOPL0382-11	Ti-MOPL0382-11
19	MOPL0382-12	Ti-MOPL0382-12
20	MOPL0382-13	Ti-MOPL0382-13
21	MOPL0382-14	Ti-MOPL0382-14
22	MOPL0382-15	Ti-MOPL0382-15

23	MOPL0382-16	Ti-MOPL0382-16
24	MOPL0382-17	Ti-MOPL0382-17
25	MOPL0382-18	Ti-MOPL0382-18
26	MOPL0382-19	Ti-MOPL0382-19
27	MOPL0382-20	Ti-MOPL0382-20
28	MOPL0382-21	Ti-MOPL0382-21
29	MOPL0382-22	Ti-MOPL0382-22
30	MOPL0382-23	Ti-MOPL0382-23
32	MOPL0382-24	Ti-MOPL0382-24
34	MOPL0382-25	Ti-MOPL0382-25
36	MOPL0382-26	Ti-MOPL0382-26
38	MOPL0382-27	Ti-MOPL0382-27
40	MOPL0382-28	Ti-MOPL0382-28
42	MOPL0382-29	Ti-MOPL0382-29
44	MOPL0382-30	Ti-MOPL0382-30
46	MOPL0382-31	Ti-MOPL0382-31
48	MOPL0382-32	Ti-MOPL0382-32
50	MOPL0382-33	Ti-MOPL0382-33


Cannulated Screw Ø 3.0mm - Short Thread

Fixation of fractures of the forearm, hand and foot, e.g.:


- Fractures and arthrodeses of the carpals and metacarpals
- Fractures of the distal radius and radial head
- Metatarsal fractures
- Drill Bit : 2.0mm



Self Drilling Tip - s.s.



Self Drilling Tip - Titanium




Code - S.S		Code - TITANIUM
Length in mm	Self Drilling	Self Drilling
8	MOPL0383-01	Ti-MOPL0383-01
9	MOPL0383-02	Ti-MOPL0383-02
10	MOPL0383-03	Ti-MOPL0383-03
11	MOPL0383-04	Ti-MOPL0383-04
12	MOPL0383-05	Ti-MOPL0383-05
13	MOPL0383-06	Ti-MOPL0383-06
14	MOPL0383-07	Ti-MOPL0383-07
15	MOPL0383-08	Ti-MOPL0383-08
16	MOPL0383-09	Ti-MOPL0383-09
17	MOPL0383-10	Ti-MOPL0383-10
18	MOPL0383-11	Ti-MOPL0383-11
19	MOPL0383-12	Ti-MOPL0383-12
20	MOPL0383-13	Ti-MOPL0383-13
21	MOPL0383-14	Ti-MOPL0383-14
22	MOPL0383-15	Ti-MOPL0383-15
23	MOPL0383-16	Ti-MOPL0383-16
24	MOPL0383-17	Ti-MOPL0383-17
25	MOPL0383-18	Ti-MOPL0383-18
26	MOPL0383-19	Ti-MOPL0383-19
27	MOPL0383-20	Ti-MOPL0383-20
28	MOPL0383-21	Ti-MOPL0383-21
29	MOPL0383-22	Ti-MOPL0383-22
30	MOPL0383-23	Ti-MOPL0383-23
32	MOPL0383-24	Ti-MOPL0383-24
34	MOPL0383-25	Ti-MOPL0383-25
36	MOPL0383-26	Ti-MOPL0383-26
38	MOPL0383-27	Ti-MOPL0383-27
40	MOPL0383-28	Ti-MOPL0383-28

42	MOPL0383-29	Ti-MOPL0383-29
44	MOPL0383-30	Ti-MOPL0383-30
46	MOPL0383-31	Ti-MOPL0383-31
48	MOPL0383-32	Ti-MOPL0383-32
50	MOPL0383-33	Ti-MOPL0383-33

Cannulated Screw Ø 3.0mm - Full Thread

→ Drill Bit : 2.0mm



Code - S.S		Code - TITANIUM
Length in mm	Self Drilling	Self Drilling
8	MOPL0384-01	Ti-MOPL0384-01
9	MOPL0384-02	Ti-MOPL0384-02
10	MOPL0384-03	Ti-MOPL0384-03
11	MOPL0384-04	Ti-MOPL0384-04
12	MOPL0384-05	Ti-MOPL0384-05
13	MOPL0384-06	Ti-MOPL0384-06
14	MOPL0384-07	Ti-MOPL0384-07
15	MOPL0384-08	Ti-MOPL0384-08
16	MOPL0384-09	Ti-MOPL0384-09


17	MOPL0384-10	Ti-MOPL0384-10
18	MOPL0384-11	Ti-MOPL0384-11
19	MOPL0384-12	Ti-MOPL0384-12
20	MOPL0384-13	Ti-MOPL0384-13
21	MOPL0384-14	Ti-MOPL0384-14
22	MOPL0384-15	Ti-MOPL0384-15
23	MOPL0384-16	Ti-MOPL0384-16
24	MOPL0384-17	Ti-MOPL0384-17
25	MOPL0384-18	Ti-MOPL0384-18
26	MOPL0384-19	Ti-MOPL0384-19
27	MOPL0384-20	Ti-MOPL0384-20
28	MOPL0384-21	Ti-MOPL0384-21
29	MOPL0384-22	Ti-MOPL0384-22
30	MOPL0384-23	Ti-MOPL0384-23
32	MOPL0384-24	Ti-MOPL0384-24
34	MOPL0384-25	Ti-MOPL0384-25
36	MOPL0384-26	Ti-MOPL0384-26
38	MOPL0384-27	Ti-MOPL0384-27
40	MOPL0384-28	Ti-MOPL0384-28
42	MOPL0384-29	Ti-MOPL0384-29
44	MOPL0384-30	Ti-MOPL0384-30
46	MOPL0384-31	Ti-MOPL0384-31
48	MOPL0384-32	Ti-MOPL0384-32
50	MOPL0384-33	Ti-MOPL0384-33

3.5mm

Cannulated Screw Ø 3.5mm Short Thread

→ Drill Bit : 2.5mm

→ Drill Bit for Gliding Hole: 3.5mm



Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
10	MOPL0385-01	Ti-MOPL0385-01
12	MOPL0385-02	Ti-MOPL0385-02
14	MOPL0385-03	Ti-MOPL0385-03
16	MOPL0385-04	Ti-MOPL0385-04
18	MOPL0385-05	Ti-MOPL0385-05
20	MOPL0385-06	Ti-MOPL0385-06
22	MOPL0385-07	Ti-MOPL0385-07
24	MOPL0385-08	Ti-MOPL0385-08
26	MOPL0385-09	Ti-MOPL0385-09
28	MOPL0385-10	Ti-MOPL0385-10
30	MOPL0385-11	Ti-MOPL0385-11
32	MOPL0385-12	Ti-MOPL0385-12

34	MOPL0385-13	Ti-MOPL0385-13
36	MOPL0385-14	Ti-MOPL0385-14
38	MOPL0385-15	Ti-MOPL0385-15
40	MOPL0385-16	Ti-MOPL0385-16
42	MOPL0385-17	Ti-MOPL0385-17
44	MOPL0385-18	Ti-MOPL0385-18
46	MOPL0385-19	Ti-MOPL0385-19
48	MOPL0385-20	Ti-MOPL0385-20
50	MOPL0385-21	Ti-MOPL0385-21

Cannulated Screw Ø 3.5mm Full Thread

→ Drill Bit : 2.5mm

→ Drill Bit for Gliding Hole: 3.5mm


Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
10	MOPL0386-01	Ti-MOPL0386-01
12	MOPL0386-02	Ti-MOPL0386-02
14	MOPL0386-03	Ti-MOPL0386-03
16	MOPL0386-04	Ti-MOPL0386-04
18	MOPL0386-05	Ti-MOPL0386-05
20	MOPL0386-06	Ti-MOPL0386-06
22	MOPL0386-07	Ti-MOPL0386-07
24	MOPL0386-08	Ti-MOPL0386-08
26	MOPL0386-09	Ti-MOPL0386-09
28	MOPL0386-10	Ti-MOPL0386-10
30	MOPL0386-11	Ti-MOPL0386-11
32	MOPL0386-12	Ti-MOPL0386-12
34	MOPL0386-13	Ti-MOPL0386-13
36	MOPL0386-14	Ti-MOPL0386-14
38	MOPL0386-15	Ti-MOPL0386-15

40	MOPL0386-16	Ti-MOPL0386-16
42	MOPL0386-17	Ti-MOPL0386-17
44	MOPL0386-18	Ti-MOPL0386-18
46	MOPL0386-19	Ti-MOPL0386-19
48	MOPL0386-20	Ti-MOPL0386-20
50	MOPL0386-21	Ti-MOPL0386-21

4.0mm

Cannulated Cancellous Screw Ø 4.0mm - Short Thread

→ Drill Bit : 2.7mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
10	MOPL0387-01	Ti-MOPL0387-01
12	MOPL0387-02	Ti-MOPL0387-02
14	MOPL0387-03	Ti-MOPL0387-03
16	MOPL0387-04	Ti-MOPL0387-04
18	MOPL0387-05	Ti-MOPL0387-05
20	MOPL0387-06	Ti-MOPL0387-06
22	MOPL0387-07	Ti-MOPL0387-07

24	MOPL0387-08	Ti-MOPL0387-08
26	MOPL0387-09	Ti-MOPL0387-09
28	MOPL0387-10	Ti-MOPL0387-10
30	MOPL0387-11	Ti-MOPL0387-11
32	MOPL0387-12	Ti-MOPL0387-12
34	MOPL0387-13	Ti-MOPL0387-13
36	MOPL0387-14	Ti-MOPL0387-14
38	MOPL0387-15	Ti-MOPL0387-15
40	MOPL0387-16	Ti-MOPL0387-16
42	MOPL0387-17	Ti-MOPL0387-17
44	MOPL0387-18	Ti-MOPL0387-18
46	MOPL0387-19	Ti-MOPL0387-19
48	MOPL0387-20	Ti-MOPL0387-20
50	MOPL0387-21	Ti-MOPL0387-21
55	MOPL0387-22	Ti-MOPL0387-22
60	MOPL0387-23	Ti-MOPL0387-23

Cannulated Cancellous Screw Ø 4.0mm - Full Thread

→ Drill Bit : 2.7mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
10	MOPL0388-01	Ti-MOPL0388-01
12	MOPL0388-02	Ti-MOPL0388-02
14	MOPL0388-03	Ti-MOPL0388-03
16	MOPL0388-04	Ti-MOPL0388-04
18	MOPL0388-05	Ti-MOPL0388-05
20	MOPL0388-06	Ti-MOPL0388-06
22	MOPL0388-07	Ti-MOPL0388-07
24	MOPL0388-08	Ti-MOPL0388-08
26	MOPL0388-09	Ti-MOPL0388-09
28	MOPL0388-10	Ti-MOPL0388-10
30	MOPL0388-11	Ti-MOPL0388-11
32	MOPL0388-12	Ti-MOPL0388-12
34	MOPL0388-13	Ti-MOPL0388-13
36	MOPL0388-14	Ti-MOPL0388-14
38	MOPL0388-15	Ti-MOPL0388-15
40	MOPL0388-16	Ti-MOPL0388-16
42	MOPL0388-17	Ti-MOPL0388-17
44	MOPL0388-18	Ti-MOPL0388-18
46	MOPL0388-19	Ti-MOPL0388-19
48	MOPL0388-20	Ti-MOPL0388-20
50	MOPL0388-21	Ti-MOPL0388-21
55	MOPL0388-22	Ti-MOPL0388-22
60	MOPL0388-23	Ti-MOPL0388-23

4.5mm

Cannulated Screw Ø 4.5mm Short Thread

→ Drill Bit : 3.2mm

→ Drill Bit for Gliding Hole: 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
20	MOPL0389-01	Ti-MOPL0389-01
22	MOPL0389-02	Ti-MOPL0389-02
24	MOPL0389-03	Ti-MOPL0389-03
26	MOPL0389-04	Ti-MOPL0389-04
28	MOPL0389-05	Ti-MOPL0389-05
30	MOPL0389-06	Ti-MOPL0389-06
32	MOPL0389-07	Ti-MOPL0389-07
34	MOPL0389-08	Ti-MOPL0389-08
36	MOPL0389-09	Ti-MOPL0389-09
38	MOPL0389-10	Ti-MOPL0389-10
40	MOPL0389-11	Ti-MOPL0389-11
42	MOPL0389-12	Ti-MOPL0389-12
44	MOPL0389-13	Ti-MOPL0389-13

46	MOPL0389-14	Ti-MOPL0389-14
48	MOPL0389-15	Ti-MOPL0389-15
50	MOPL0389-16	Ti-MOPL0389-16
52	MOPL0389-17	Ti-MOPL0389-17
54	MOPL0389-18	Ti-MOPL0389-18
56	MOPL0389-19	Ti-MOPL0389-19
58	MOPL0389-20	Ti-MOPL0389-20
60	MOPL0389-21	Ti-MOPL0389-21
62	MOPL0389-22	Ti-MOPL0389-22
64	MOPL0389-23	Ti-MOPL0389-23
66	MOPL0389-24	Ti-MOPL0389-24
68	MOPL0389-25	Ti-MOPL0389-25
70	MOPL0389-26	Ti-MOPL0389-26
72	MOPL0389-27	Ti-MOPL0389-27
74	MOPL0389-28	Ti-MOPL0389-28

Cannulated Screw Ø 4.5mm - Full Thread

→ Drill Bit : 3.2mm

→ Drill Bit for Gliding Hole: 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
20	MOPL0390-01	Ti-MOPL0390-01
22	MOPL0390-02	Ti-MOPL0390-02
24	MOPL0390-03	Ti-MOPL0390-03
26	MOPL0390-04	Ti-MOPL0390-04
28	MOPL0390-05	Ti-MOPL0390-05
30	MOPL0390-06	Ti-MOPL0390-06
32	MOPL0390-07	Ti-MOPL0390-07
34	MOPL0390-08	Ti-MOPL0390-08
36	MOPL0390-09	Ti-MOPL0390-09
38	MOPL0390-10	Ti-MOPL0390-10
40	MOPL0390-11	Ti-MOPL0390-11
42	MOPL0390-12	Ti-MOPL0390-12
44	MOPL0390-13	Ti-MOPL0390-13
46	MOPL0390-14	Ti-MOPL0390-14
48	MOPL0390-15	Ti-MOPL0390-15
50	MOPL0390-16	Ti-MOPL0390-16
52	MOPL0390-17	Ti-MOPL0390-17
54	MOPL0390-18	Ti-MOPL0390-18
56	MOPL0390-19	Ti-MOPL0390-19
58	MOPL0390-20	Ti-MOPL0390-20
60	MOPL0390-21	Ti-MOPL0390-21
62	MOPL0390-22	Ti-MOPL0390-22
64	MOPL0390-23	Ti-MOPL0390-23
66	MOPL0390-24	Ti-MOPL0390-24
68	MOPL0390-25	Ti-MOPL0390-25
70	MOPL0390-26	Ti-MOPL0390-26
72	MOPL0390-27	Ti-MOPL0390-27
74	MOPL0390-28	Ti-MOPL0390-28

6.5mm

Cancellous Cannulated Screw Ø 6.5mm 16mm Thread

→ Drill Bit : 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
30	MOPL0391-01	Ti-MOPL0391-01
35	MOPL0391-02	Ti-MOPL0391-02
40	MOPL0391-03	Ti-MOPL0391-03
45	MOPL0391-04	Ti-MOPL0391-04
50	MOPL0391-05	Ti-MOPL0391-05
55	MOPL0391-06	Ti-MOPL0391-06
60	MOPL0391-07	Ti-MOPL0391-07
65	MOPL0391-08	Ti-MOPL0391-08
70	MOPL0391-09	Ti-MOPL0391-09
75	MOPL0391-10	Ti-MOPL0391-10
80	MOPL0391-11	Ti-MOPL0391-11
85	MOPL0391-12	Ti-MOPL0391-12
90	MOPL0391-13	Ti-MOPL0391-13
95	MOPL0391-14	Ti-MOPL0391-14
100	MOPL0391-15	Ti-MOPL0391-15

105	MOPL0391-16	Ti-MOPL0391-16
110	MOPL0391-17	Ti-MOPL0391-17
115	MOPL0391-18	Ti-MOPL0391-18
120	MOPL0391-19	Ti-MOPL0391-19
125	MOPL0391-20	Ti-MOPL0391-20
130	MOPL0391-21	Ti-MOPL0391-21
135	MOPL0391-22	Ti-MOPL0391-22
140	MOPL0391-23	Ti-MOPL0391-23
145	MOPL0391-24	Ti-MOPL0391-24
150	MOPL0391-25	Ti-MOPL0391-25

Cancellous Cannulated Screw Ø 6.5mm 32mm Thread

→ Drill Bit : 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
40	MOPL0392-01	Ti-MOPL0392-01
45	MOPL0392-02	Ti-MOPL0392-02
50	MOPL0392-03	Ti-MOPL0392-03
55	MOPL0392-04	Ti-MOPL0392-04

60	MOPL0392-05	Ti-MOPL0392-05
65	MOPL0392-06	Ti-MOPL0392-06
70	MOPL0392-07	Ti-MOPL0392-07
75	MOPL0392-08	Ti-MOPL0392-08
80	MOPL0392-09	Ti-MOPL0392-09
85	MOPL0392-10	Ti-MOPL0392-10
90	MOPL0392-11	Ti-MOPL0392-11
95	MOPL0392-12	Ti-MOPL0392-12
100	MOPL0392-13	Ti-MOPL0392-13
105	MOPL0392-14	Ti-MOPL0392-14
110	MOPL0392-15	Ti-MOPL0392-15
115	MOPL0392-16	Ti-MOPL0392-16
120	MOPL0392-17	Ti-MOPL0392-17
125	MOPL0392-18	Ti-MOPL0392-18
130	MOPL0392-19	Ti-MOPL0392-19
135	MOPL0392-20	Ti-MOPL0392-20
140	MOPL0392-21	Ti-MOPL0392-21
145	MOPL0392-22	Ti-MOPL0392-22
150	MOPL0392-23	Ti-MOPL0392-23

Cancellous Cannulated Screw Ø 6.5mm Full Thread

→ Drill Bit : 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
30	MOPL0393-01	Ti-MOPL0393-01
35	MOPL0393-02	Ti-MOPL0393-02
40	MOPL0393-03	Ti-MOPL0393-03
45	MOPL0393-04	Ti-MOPL0393-04
50	MOPL0393-05	Ti-MOPL0393-05
55	MOPL0393-06	Ti-MOPL0393-06
60	MOPL0393-07	Ti-MOPL0393-07
65	MOPL0393-08	Ti-MOPL0393-08
70	MOPL0393-09	Ti-MOPL0393-09
75	MOPL0393-10	Ti-MOPL0393-10
80	MOPL0393-11	Ti-MOPL0393-11
85	MOPL0393-12	Ti-MOPL0393-12
90	MOPL0393-13	Ti-MOPL0393-13
95	MOPL0393-14	Ti-MOPL0393-14
100	MOPL0393-15	Ti-MOPL0393-15
105	MOPL0393-16	Ti-MOPL0393-16

110	MOPL0393-17	Ti-MOPL0393-17
115	MOPL0393-18	Ti-MOPL0393-18
120	MOPL0393-19	Ti-MOPL0393-19
125	MOPL0393-20	Ti-MOPL0393-20
130	MOPL0393-21	Ti-MOPL0393-21
135	MOPL0393-22	Ti-MOPL0393-22
140	MOPL0393-23	Ti-MOPL0393-23
145	MOPL0393-24	Ti-MOPL0393-24
150	MOPL0393-25	Ti-MOPL0393-25

7.0mm

Cancellous Cannulated Screw Ø 7.0mm, 16mm Thread

→ Drill Bit : 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
50	MOPL0394-01	Ti-MOPL0394-01
55	MOPL0394-02	Ti-MOPL0394-02
60	MOPL0394-03	Ti-MOPL0394-03

65	MOPL0394-04	Ti-MOPL0394-04
70	MOPL0394-05	Ti-MOPL0394-05
75	MOPL0394-06	Ti-MOPL0394-06
80	MOPL0394-07	Ti-MOPL0394-07
85	MOPL0394-08	Ti-MOPL0394-08
90	MOPL0394-09	Ti-MOPL0394-09
95	MOPL0394-10	Ti-MOPL0394-10
100	MOPL0394-11	Ti-MOPL0394-11
105	MOPL0394-12	Ti-MOPL0394-12
110	MOPL0394-13	Ti-MOPL0394-13
115	MOPL0394-14	Ti-MOPL0394-14
120	MOPL0394-15	Ti-MOPL0394-15

Cancellous Cannulated Screw Ø 7.0mm, 32mm Thread

→ Drill Bit : 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
50	MOPL0395-01	Ti-MOPL0395-01

55	MOPL0395-02	Ti-MOPL0395-02
60	MOPL0395-03	Ti-MOPL0395-03
65	MOPL0395-04	Ti-MOPL0395-04
70	MOPL0395-05	Ti-MOPL0395-05
75	MOPL0395-06	Ti-MOPL0395-06
80	MOPL0395-07	Ti-MOPL0395-07
85	MOPL0395-08	Ti-MOPL0395-08
90	MOPL0395-09	Ti-MOPL0395-09
95	MOPL0395-10	Ti-MOPL0395-10
100	MOPL0395-11	Ti-MOPL0395-11
105	MOPL0395-12	Ti-MOPL0395-12
110	MOPL0395-13	Ti-MOPL0395-13
115	MOPL0395-14	Ti-MOPL0395-14
120	MOPL0395-15	Ti-MOPL0395-15

Cancellous Cannulated Screw Ø 7.0mm, Full Thread

→ Drill Bit : 4.5mm




Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
50	MOPL0396-01	Ti-MOPL0396-01
55	MOPL0396-02	Ti-MOPL0396-02
60	MOPL0396-03	Ti-MOPL0396-03
65	MOPL0396-04	Ti-MOPL0396-04
70	MOPL0396-05	Ti-MOPL0396-05
75	MOPL0396-06	Ti-MOPL0396-06
80	MOPL0396-07	Ti-MOPL0396-07
85	MOPL0396-08	Ti-MOPL0396-08
90	MOPL0396-09	Ti-MOPL0396-09
95	MOPL0396-10	Ti-MOPL0396-10
100	MOPL0396-11	Ti-MOPL0396-11
105	MOPL0396-12	Ti-MOPL0396-12
110	MOPL0396-13	Ti-MOPL0396-13
115	MOPL0396-14	Ti-MOPL0396-14
120	MOPL0396-15	Ti-MOPL0396-15

7.3mm

Cancellous Cannulated Screw Ø 7.3mm - 16mm Thread


→ Drill Bit : 5.0mm



Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
50	MOPL0397-01	Ti-MOPL0397-01
55	MOPL0397-02	Ti-MOPL0397-02
60	MOPL0397-03	Ti-MOPL0397-03
65	MOPL0397-04	Ti-MOPL0397-04
70	MOPL0397-05	Ti-MOPL0397-05
75	MOPL0397-06	Ti-MOPL0397-06
80	MOPL0397-07	Ti-MOPL0397-07
85	MOPL0397-08	Ti-MOPL0397-08
90	MOPL0397-09	Ti-MOPL0397-09
95	MOPL0397-10	Ti-MOPL0397-10
100	MOPL0397-11	Ti-MOPL0397-11
105	MOPL0397-12	Ti-MOPL0397-12
110	MOPL0397-13	Ti-MOPL0397-13
115	MOPL0397-14	Ti-MOPL0397-14
120	MOPL0397-15	Ti-MOPL0397-15

Cancellous Cannulated Screw Ø 7.3mm - 32mm Thread


→ Drill Bit : 5.0mm



Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
50	MOPL0398-01	Ti-MOPL0398-01
55	MOPL0398-02	Ti-MOPL0398-02
60	MOPL0398-03	Ti-MOPL0398-03
65	MOPL0398-04	Ti-MOPL0398-04
70	MOPL0398-05	Ti-MOPL0398-05
75	MOPL0398-06	Ti-MOPL0398-06
80	MOPL0398-07	Ti-MOPL0398-07
85	MOPL0398-08	Ti-MOPL0398-08
90	MOPL0398-09	Ti-MOPL0398-09
95	MOPL0398-10	Ti-MOPL0398-10
100	MOPL0398-11	Ti-MOPL0398-11
105	MOPL0398-12	Ti-MOPL0398-12
110	MOPL0398-13	Ti-MOPL0398-13
115	MOPL0398-14	Ti-MOPL0398-14
120	MOPL0398-15	Ti-MOPL0398-15

Cancellous Cannulated Screw Ø 7.3mm - Full Thread

→ Drill Bit : 5.0mm



Stainless Steel		Titanium
Length in mm	Self Drilling	Self Drilling
50	MOPL0399-01	Ti-MOPL0399-01
55	MOPL0399-02	Ti-MOPL0399-02
60	MOPL0399-03	Ti-MOPL0399-03
65	MOPL0399-04	Ti-MOPL0399-04
70	MOPL0399-05	Ti-MOPL0399-05
75	MOPL0399-06	Ti-MOPL0399-06
80	MOPL0399-07	Ti-MOPL0399-07
85	MOPL0399-08	Ti-MOPL0399-08
90	MOPL0399-09	Ti-MOPL0399-09
95	MOPL0399-10	Ti-MOPL0399-10
100	MOPL0399-11	Ti-MOPL0399-11
105	MOPL0399-12	Ti-MOPL0399-12
110	MOPL0399-13	Ti-MOPL0399-13
115	MOPL0399-14	Ti-MOPL0399-14
120	MOPL0399-15	Ti-MOPL0399-15



MADIHA ORTHO

Solution for humanity



support@madahaortho.com

info@madahaortho.com

sales@madahaortho.com

Contact: +918448115377
+917003006585

WhatsApp: +918697716553
+917602449449

www.madahaortho.com

Office: 542,4th floor, G.T. Road, Howrah-71101
(West Bengal, India)